

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# CCNG1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag1157

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MI EVLT TTD SQ KLL HQ L NALLE QES RC QPK VC GL RL IE  
SA HDNGL RMT AR LR DF EV KD LLS LT Q FF GF DT ET FS LA  
VN LLD RF LSK MK VQ PK HL GC VGL SC FY LA VKS IEE ERN  
V PLAT DL IRIS QY RF TV SD LMR MEK I VLE KVC WK VK AT  
T AF Q FL QLY YS LL QEN LPL ER RN SIN FER LEA QLK ACH  
CR IIF SK AK PS VL ALS IIA LE IQA QK CVEL TE GIE CL QKH  
SK INGR DLT FWQ ELV SK CL TE YSS NE CS KPN VQ K LKW  
IV SGR TAR QL KHS YYR ITH LPT IPE MVP  
(1-295 aa encoded by BC007093)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.