

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human Synaptophysin; SYP protein



Basic Information

Catalog Number:

Ag11781

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ELQLSVDCA~~NKTES~~DL~~SIEVE~~F~~EY~~P~~FRL~~HQ~~VY~~F~~DAP~~TC
RGGTTK~~VFL~~VG~~DYSS~~AE~~FFVT~~VAV~~FAFL~~YS~~MG~~ALAT
YIFLQ~~NKYRE~~NN~~KGP~~ML~~DFLA~~TAV~~FAFM~~WL~~VSS~~SAW
AKGLSD~~VKMAT~~DPENII~~KEMP~~VCRQT~~GN~~TC~~KEL~~RDLV
TSG~~LNTSV~~V~~FGFL~~N~~LVL~~W~~VGNL~~W~~FWF~~K~~ETG~~WAAPFL
RAP~~PGA~~PEKQ~~PAPG~~DAY~~GDAG~~Y~~GQGP~~GGY~~GPQ~~DSY
GPQGGY~~QPDY~~GQ~~PAGS~~GG~~SGYGP~~Q~~DYGQ~~Q~~QGY~~GP
QGAP~~TSFS~~NQM

(50-313 aa encoded by BC064550)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.