

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human PLEKHA9 protein



## Basic Information

Catalog Number:

Ag11854

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSELRLCCDLLVQQVDKTEVTTTGVSNSEEGIDVGT  
LLKSTCNTFLKTEECMQIANAAFTSELYHTPPGSPQ  
LAMLKSSKMKHPHPIHNSLERQTELSTCENGSLNMEIN  
GEEIILMNKNSLYLKSIEDCSISSEENTDDNITVQGE  
IMKEDRMENLKNHDNNLSQSGSDSSCSPECLWEEGK  
EVIPTFFSTMNTSFSIDIELLEDGIPTEAFASCCAVVP  
VLDKLGPTVFAPVKMDLVENIKKVNQKYITNKEEFTTL  
QKIVLHEVEADVAQVRNSATEALLWLKRGKFLKGF  
TEVKNGEKDIQTALNNAYGKTLRQHHGWVVRGVFAL  
ALRATPSYEDFVAALTVKEGDHRKEAFSIGMQRDLSL  
YLPAMKKQMAILDAL

(1-391 aa encoded by BC063575)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.