

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human CD146, MCAM protein



## Basic Information

### Catalog Number:

Ag11855

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Peptide Sequence:

GEAEQPAPELVEVEVGSTALLKCGLSQSQGNLSHVD  
WFSVHKEKRTLIFRVRQGQGGQSEPGYEQRSLQDR  
GATLALTQVTPQDERIFLCQGKRPRSQYRIQLRVYK  
APEEPNIQVNPLGIPVNSKEPEEVATCVGRNGYPIQ  
VIWYKNGRPLKEEKNRVHIQSSQTVESGLYTLQSILK  
AQLVKEDKDAQFYCELNYRLPSGNHMKESREVTVPV  
FYPTKVVLEVPVGMLKEGDRVEIRCLADGNPPPHF  
SISKQNPSTREAEETTNDNGVLVLEPARKEHSGRYEC  
QGLDLDTMISLLSEPQELLVNYVSDVRVSP  
(26-350 aa encoded by BC056418)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.