

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human THNSL2 protein



Basic Information

Catalog Number:

Ag12020

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6^oHis.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MWVYSTRGVAPRVNFEALFSGYAPDGGLFMPEELP
QLDRETLQCQWSTLSYPGLVKELCALFIGSELLPKDELN
DLIDRAFSRFRHREVVHLSRLRNLNVLELWHGVTYAF
KDLSLCTTQFLQYFLEKREKHVTVVVGTSGDTGSAA
IESVQGAKNMDIIVLLPKGHCTKIQELQMTTVLKQNV
HVFVGEVGNDELDEPIKTVFADVAFVKKHNLMSLNSI
NWSRVLVQMAHHFFAYFQCTPSLDTHPLPLVEVVVP
TGAAGNLAAGYIAQKIGLPIRLVAVVNRNDIIHRTVQ
QGDFSLSEAVKSTLASAMDIQVPYNMEGVFWLLSGS
DSQVTRALMEQFERTQSVNLPKELHSHKWERQDYKEM
AVMECDGCCVELCLGNCGPRRGSVTDIPGTAMVRG
GG

(1-405aa encoded by BC064423)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.