

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ERRa protein



Basic Information

Catalog Number:
Ag12042

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MSSQVVGIEPLYIKAEPASPDSPKGSSETETEPVALA
PGPAPTRCLPGHKEEEDGEGAGPGEQGGGKLVLSLP
KRLCLVCGDVASGYHYGVASCEACKAFFKRTIQGSIE
YSCPASNECEITKRRRKACQACRFTKLRVGMKEGV
RLDRVRGGRQKYKRRPEVDPLPFGPFPAAGPLAVAGG
PRKTAAPVNALVSHLLVVEPEKLYAMPDPAAGPDGHL
AVATLCDLFDREIVVTISWAKSIPGFSSLSLSDQMSVL
QSVWMEVLVLGVAQRSLPLQDELAF AEDLVLDEEGA
RAAGLGEALGALLQLVRRQLRLEREYVLLKALALA
NSDSVHIEDAEAVEQLREALHEALLEYEAGRAGPGGG
AERRRAGRLLLTPLLRQTAGKVLAHFYGVKLEGKVP
MKLFLEMLEAMMD
(1-423aa encoded by BC092470)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W:
(toll free in USA), or 1(312) 455-8498 ptglab.com
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.