

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ADAL protein



Basic Information

Catalog Number:

Ag12114

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MI E A E E Q Q P C K T D F Y S E L P K V E L H A H L N G S I S S H T M K
K L J A Q K P D L K I H D Q M T V I D K G K K R T L E E C F Q M F Q T I H
Q L T S S P E D I L M V T K D V I K E F A D D G V K Y L E R S T P R R E N
A T G M T K K T Y V E S I L E G I K Q S K Q E N L I D V R Y L I A V D R R
G G P L V A K E T V K L A E E F F L S T E G T V L G L D L S G D P T V G Q
A K D F L E P L L E A K I A G L K L A L H L S E I P N Q K K E T Q I L L D L L P
D R I G H G T F L N S G E G G S L D L V D F V R Q H R I P L G K A W S F R
S S R

(1-267 aa encoded by BC075857)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.