

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# KLHL32 Fusion Protein



## Basic Information

**Catalog Number:**

Ag12196

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MPSERCLSIQEMLTGQRLCHSESHNDSVLAALNQQR  
SDGILCDITLIAEEQKFHAHKAVLAACSDYFRAMFSLC  
MVESGADEVNLHGVTSLGLKQALEFAYTGQILLEPGV  
IQDVLAAAGSHLQLELLNLCSHYLIQELNSFNLDLYRL  
ADLFNLTLEKAVIDFLVKHLSSELLKSRPEEVLTPYCLL  
QEVLKSDRLTSLSEEQIWQLAVRWLEHNCYQYMDE  
LLQYIRFGLMDVDTLHTVALSHPLVQASETATLVNE  
ALEYHQSIYAQPVVWQTRRTKPRFQSDTLYIIGGKKRE  
VCKVKELRYFNPVDQENALIAAIANWSELAPMPVGRS  
HHCVAVMGDFLFVAGGEVEHASGRTCAVRTACRYD  
PRNSWAEIAPMKNCREHFVLGAMEEYLYAVGGRNE  
LRQVLPTVERYCPKKNWTFVQSFDRSLSCHAGYVA  
DGLLWISGRTYLMLDLSKHTFIVVYI  
(1-470 aa encoded by BC071627)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.