

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human INPP5F protein



Basic Information

Catalog Number:

Ag12323

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

VKQTKSNVNI~~GNLRKLG~~NFTKPEMKV~~NFLKPNL~~KVNL
WKS~~DSSLET~~MENTGVMDKVQAESD~~GDMS~~SDNDSYH
SDEFLTNSK~~SDEDRQLANS~~LESVGPIDYVLPSCGIIAS
APRLGSR~~SQSL~~SST~~DSSV~~HAPSEITVAHGSLGK~~GQE~~
SPLK~~KSPS~~AGDVHILTGF~~AKPMDI~~YCHRFVQDAQNKV
THLSE~~TRSV~~SQQASQERNQMTNQV~~SD~~ETQSESTEQ
TPSR~~PSQL~~DVSLSATG~~PQFL~~SVEPAHSVA~~SQK~~TP~~SA~~
SSM~~LE~~LETGLHVT~~PSPE~~SSSR~~AVSP~~FAKIRSSMVQV
ASITQAGLTHGINFAVSKVQ~~KSPPE~~PIINQVQ~~NEL~~
KKMFIQCQTRIIQI

(C-term-345aa encoded by BC111493)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.