

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human NUP155 protein



Basic Information

Catalog Number:

Ag12388

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

NNPTVSGMSDMDYPLQGPGLLSVPNLPEISSIRRVPLP
PELVEQFGHMQCNCMMGVFPPISRAWLTIDSDIFM
WNYEDGGDLAYFDGLSETILAVGLVKPKAGIFQPHVR
HLLVLATPVDIVILGLSYANLQTGSGVLNDSLSGGMQL
LPDPLYSLPTDNTYLLTITSTDNGRIFLAGKDGCLYEVA
YQAEAGWFSQRCKINHSSLSFLVPSLLQFTFSED
DPILQIAIDNSRNILYRSEKGVIVYDLGQDGGMSR
VASVSQNAIVSAAGNIARTIDRSVFKPIVQIAVIENSE
SLDCQLLAVTHAGVRLYFSTCPFRQPLARPNTLTVH
VRLPPGFSASSTVEKPSKVHRLYSKGILLMAASENED
NDILWCNVHDTFFPQKPMMETQMTAGVD
(53-455aa encoded by BC039257)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.