

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

AURKB Fusion Protein



Basic Information

Catalog Number:
Ag1246

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
QIQLFAGGRESSALDPSSPPFSLRMAQKENSYPWP
YGRQTAPSGLSTLPQRVLRKEPVTPSALVLMRSNVQ
PTAAPGQKVMENSSGTPDILTRRHFTIDDFEIGRPLGK
GKFGNVYLAREKKSHFIVALKVLFKSQIEKEGVEHQLR
REIEIQAHLHHPNILLRNYFYDRRRYLILEYAPRGELYK
ELQKSCTFDEQRTATIMEELADALMYCHGKVIHRDI
KPENLLGLKELKIADFGWSVHAPSLRRKTMCGTLDY
LPPEMIEGRMHNEKVDLWICIVLCYELLVGNPPFESA
SHNETYRRIKVDLKFPAVPMGAQDLISKLLRHPSE
RLPLAQVSAHPWVRANSRRVLPSPALQSV
(1-371 aa encoded by BC013300)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.