

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human FASTKD3 protein



## Basic Information

**Catalog Number:**

Ag12959

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

PLI I K L G K Y V V R H V P H F T N E E L R R V L E A F I Y F G H H D T F F  
T K A L E H R V A A V C L T D P E V V C R V M E Y C S R E L I L S K P I L  
N A V A E T F V C Q T E K F S P R Q I S A L M E P F G K L N Y L P P N A S  
A L F R K L E N V L F T H F N Y P P K S L L K L L H S C S L N E C H P V N F  
L A K I F K P L F L Q R L Q G K E S H L D T L S R A Q L T Q L F L A S V L E C  
P F Y K G P K L L P K Y Q V K S F L T P C C S L E T P V D S Q L Y R Y V K I G  
L T N L L G A R L Y F A P K V L T P Y C Y T I D V E I K L D E E G F V L P S T  
A N E D I H K R I A L C I D G P K R F C S N S K H L L G K E A I K Q R H L Q L  
L G Y Q V V Q I P Y H E I G M L K S R R E L V E Y L Q R K L F S Q N T V H  
W L Q E

(C-term-350aa encoded by BC001295)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.