

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human NUP155 protein



## Basic Information

Catalog Number:

Ag12974

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, T-HIS, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

HGPEEDIVGLQAFQERLNSYKCIDTLQELVNQSKA  
APQSPSPKPKPGPPVLSDDPNMLSNEEAGHHFEQML  
KLSQRSKDELFSIALYNWLIQVDLADKLLQVASPFLEP  
HLVRMAKVDQNRVRYMDLLWRYYEKNRSFSNAARVL  
SRLADMHSTEISLQQRLEYIARAILSAKSSTAISIAAD  
GEFLHELEEKMEVARIQLQIQETLQRQYSHHSSVQDA  
VSQLDSELMIDITKLYGEFADPFKLAECKLAIHCAGYS  
DPILVQTLWQDIIEKELSDSVTLSSSDRMHALSLKIVLL  
GKIYAGTPRFPLDFIVQFLEQQVCTLNWDVGFVIQT  
MNEIGVPLRLLLEVYDQLFKSRDPFWNRMKKPLHLLDC  
IHVLLIRYVENPSQVLCERRRFTNLCLDAVCGYLVEL  
(952-1391aa encoded by BC039257)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.