

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human RUNX3 protein



Basic Information

Catalog Number:
Ag13066

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MASNSIFDSFPTYSPFIRDPSTSRFFPPSPAFCGG
GGGKMGENSGALSAQA AVGPGGARPEVRSMVDVL
ADHAGELVRTDSPNFLCSVLP SHWRCNKTLPVAFKV
VALGDVDPDGT VVTVMAGNDENYSAELRNASAVMKN
QVARFNDLRFVGRSGRGKSFLLTITVFTNPTQVATYH
RAIKVTVDGPREPRRHQKLEDQTKPFDRFGDLERLR
MRVTPSTPSRGSLSLTS HFSSQPQTPIQGTSELNPF5
DPRQFDRSFTLPTLTESRF PDPRMHYPGAMSAAFPY
SATPSGTSISSLSVAGMPATSRFHHTYLPPYPGAPQ
NQSGPFQANPSYHLYYGTSSGSYQFSMVAGSSSGG
DRSPTRMLASCTSSAASVAAGNLMNPSLGGQSDGVE
(1-429 aa encoded by BC013362)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.