

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant VP7 protein



Basic Information

Catalog Number:

Ag13422

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MDTIAARALTMACATLQEARIVLEANVMEILGIAIN
RYNGLTLRGVTMRPTSLAQRNEMFFMCLDMMLSAAG
INVGPISPDYTQHMATIGVLATPEIPFTTEAANEIARV
TGETSTWGPARGPYGFFLETEETFQPGRWFMRRAQ
AVTAVVCGPDMIQVSLNAGARGDVQQIFQGRNDP
MMIYLVWRRIFENFAMAQGNSQQTQAGVTVSVGGV
DMRAGRIIAWDGQAALHVNPTQQNAMVQIQVVF
YISMDKTLNQYPALTAEIFNVYSFRDHTWHGLRTAILN
RTTLPNMLPPIFPNDRDSILTLLLLSTLADVYTVLRPEF
AIHGVNPMPLTRAIARAAYV
(1-349aa encoded by NC_006022)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.