

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human CKAP5 protein



## Basic Information

Catalog Number:

Ag13464

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

NI L S A L L V L L Q D S L L A T A S S P K F S E L V M K L W R M V R L L P  
D T I N S I N L D R I L L D I H I F M K V F P K E K L Q C K S E F P I R T L K T  
L L H T L C K L K G P K I L D H L T M I D N K N E S E L E A H L C R M M K H  
S M D Q T G S K S D K E T E K G A S R I D E K S S K A K V N D F L A E I F K  
K I G S K E N T K E G L A E L Y E Y K K K Y S D A D I E P F L K N S S Q F F Q  
S Y V E R G L R V I E M E R E G I G R I S T S T G I S P Q M E V T C V P T P  
T S T V S S I G N T N G E E V G P S V Y L E R L K I L R Q R C G L D N T K Q  
D D R P P L T S L L S K P A V P T V A S S T D M L H S K L S Q L R E S R E Q  
H Q H S D L D S N Q T H S S G T V T S S S T A N I D D L K K R L E R I K S  
S R K

(1683-2032aa encoded by BC035554)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.