

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ADH1C Fusion Protein



## Basic Information

### Catalog Number:

Ag13528

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Validated Application:

Blocking peptide

### Peptide Sequence:

MSTAGKVIKCKAAVLWELKPFSEIEVEVAPPKAHEV  
RIKMVAAGICRSDEHVVSGLNLTPLPVILGHEAAGIVE  
SVGEGVTTVKPGDKVPLFTPQCGKCRICKNPESNYC  
LKNDLGNPRGTLQDGTTRFTCSGKPIHHFVGVSTFSQ  
YTVVDENAVAKIDAASPLEKVCLIGCGFSTGYGSAVK  
VAKVTPGSTCAVFGGLGGVGLSVVMGCKAAGAARIA  
VDINKDKFAKAKELGATECINPDYKKPIQEVLKEMT  
DGGVDFSEVIGQLDTMMASLLCCEACGTSVIVGV  
PPDSQNLNINPMLLTGRTWKGAIFGGFKSKESVPKLV  
ADFMKKFSLDALITNVLPEKINEGFDLLRSGKSIRTV  
LT

(1-374 aa encoded by BC062476)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.