

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human LINGO1 protein



## Basic Information

### Catalog Number:

Ag13597

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Peptide Sequence:

SGSATGCPPRCECSAQDRAVLCHRRFVAVPEGIPT  
TRLLDLGKNRIKTLNQDEFASFPHLEELNENIVSAVE  
PGAFNNLFNLRITGLRSNRLKLIPLGVFTGLSNLTKLDIS  
ENKIVILLDYMFQDLYNLKSLEVGDNDLVYISHRAFSG  
NSLEQLTLEKCNLTSPTEALSHLHGLVLRRLNINAI  
RDYSFKRLYRLKVLEISHWPYLDTPNCLYGLNLTSL  
ITHCNLTAVPYLAVRHLVYLRFLNLSYNPISTIEGSM  
ELLRLQEIQLVGGQLAVVEPYAFRGLNYLRVLNVSG  
QLTTLE

(36-351 aa encoded by BC011057)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.