

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human KCNS1 protein



Basic Information

Catalog Number:

Ag13608

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MLMLLVRGTHYENLRISKVVLPTPLGGRSTETVSEFPG
PDTGIRWRRSDEALRVNVGGVRRQLSARALARFPGTR
LGR LQAAA S E E Q A R R L C D D Y D E A A R E F Y F D R H P G F F L
S L L H F Y R T G H L H V L D E L C V F A F G Q E A D Y W G L G E N A L A
A C C R A R Y L E R R L T Q P H A W D E D S D T P S S V D P C P D E I S D
V Q R E L A R Y G A A R C G R L R R R L W L T M E N P G Y S L P S K
(1-220aa encoded by BC075033)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.