

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# INSR Fusion Protein



## Basic Information

Catalog Number:

Ag13922

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>His</sup>.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

RKRQPDGPLGPLYASSNPEYLSASDVFPSCSVYVPDEW  
EVSREKITLLRELQGSFGMVYEGNARDIIKGEAETRV  
AVKTVNESASLRERIEFLNEASVMKGFCHHVVRLG  
VVSKGQPTLVVMEMLAHGDLKSYLRSLRPEAENNPGR  
PPPTLQEMIQMAAEIADGMAYLNAKFFVHRDLAARN  
CMVAHDFTVKGDFGMTRDIYETDYRKGKGLLPVR  
WMAPESLKDGVFTTSSDMWSFGVVLWEITSLAEQPY  
QGLSNEQVLKFMVMDGGYLDQPDNCPERVTDLMRMC  
WQFNPKMRPTFLEIVNLLKDDLHPSFPEVSFFHSEENK  
APESELEMEFEDMENVPLDRSSHQCREEAGGRDGG  
SSLGFKRSYEEHIPYTHMNGGKKNRILTLPRSNPS  
(968-1370 aa encoded by BC117172)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.