

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

OSMR Fusion Protein



Basic Information

Catalog Number:
Ag1426

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MALFAVFQTTFFLTLRLTYQSEVLAERLPLTPVSLKV
STNSTRQSLHLQWTVHNLPHYQELKMFQIQISRIET
SNVIWVGNYSTTVKWNQVLHWSWESELPLECATHF
VRIKSLVDDAKFPEPNFWSNWSSWEEVSVQDSTGQ
DILFVFPKDKLVEEGTNVTICYVSRNIQNNVSCYLEGK
QIHGEQLDPHVTAFNLNSVPPFIRNKGNTNIYCEASQGN
VSEGMKGIVLFVSKVLEEPKDFSCETEDFKLHCTWD
PGTDTALGWSKQPSQSYTLFESFSGEKKLCTHKNWC
NWQITQDSQETYNFTLIAENYLKRSVNILFNLTHRGE
TRVVTARHG
(1-342 aa encoded by BC010943)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.