

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human FES protein



## Basic Information

**Catalog Number:**

Ag14641

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

VAELLVHSGDFLVRESQGKQEVVLSVLWDGLPRHFII  
QSLDNLYRLEGEGFPSIPLLDHLLSTQQPLTKKSGVVL  
HRAVPKDKWWLNHEDLVLGEQIGRGNFGEVFSGRLLR  
ADNTLVAVKSCRETLPDLKAKFLQEARILKQYSHPNL  
VRLIGVCTQKQPIYIVMELVQGGDFLTLRTEGARLRV  
KTLQMVGDAAAGMEYLESKCCIHRDLAARNCLVTEK  
NVLKISDFGMSREEADGVYAAAGGLRQVPVKWTAPE  
ALNYGRYSSESDVWSFGILLWETFSLGASYPNLSNQ  
QTREFVEKGGRLPCPELCPDAVFRLMEQCWAYEPGQ  
RPSFSTIYQELQSIKRRHR

(C-term-351aa encoded by BC035357)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.