

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human IL-1F5 protein



## Basic Information

Catalog Number:

Ag14658

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MVLSGALCFRMKDSALKVLYLHNNQLLAGGLHAGKVI  
KGEEISVVPNRWLDASLSPVILGVQGGSQCLSCGVG  
QEPTLTLEPVNIMELYLGAKESKSFTFYRRDMGLTSSFE  
SAAYPGWFLCTVPEADQPVRTQLPENGGWNPITD  
FYFQQCD

(1-155 aa encoded by BC024747)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.