

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TNK2 Fusion Protein



## Basic Information

Catalog Number:

Ag1472

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MQPEEGTGWLLELLSEVQLQQYFLRLRDDLNVTRLSH  
FEYVKNEDLEKIGMGRPGQRRLWEAVKRRKALCKRKS  
WMSKVFSGKRLEAEFPPHHSQSTFRKTSPAPGGPAG  
EGPLQSLTCLIGEKDLRLEKLGDSFVVVRRGEWDAP  
SGKTVSVAVKCLKPDVLSQPEAMDDFIREVNAMHSL  
DHRNLIIRLYGVVLTTPMKMVTTELAPLGSLLDRLRKHQG  
HFLLGTLRYAVQVAEGMGYLESKRFIHRDLAARNLLL  
ATRDLVKIGDFGLMRALPQNDDHYVMQEHKVPFA  
WCAPE

(N-term-300aa encoded by BC008884)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.