

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human BRAF protein



Basic Information

Catalog Number:

Ag15014

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

NEVGVLKTRHVNILLFMGYSTKPQLAIVTQWCEGSS
LYHHLHIIETKFEMIKLIDARQTAQGM DY L H A K S I I H R
DLKSNNIFLHEDLTVKIGDFGLATVKSRWSGSHQFEQ
LSGSILWMAPEVIRMQDKNPYSFQSDVYAFGIVLYEL
MTGQLPYSNINNRDQIIFMVGRGYLSPDLKVRSNCP
KAMKRLMAECLKKKRDERPLFPQILASIELARSLPKIHR
SASEPSLNRAGFQTEDFSLYACASPKTPIQAGGYGAF
PVH

(500-766 aa encoded by BC101757)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.