

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human TIGD7 protein



## Basic Information

Catalog Number:

Ag15597

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

DI~~C~~LP~~G~~KK~~I~~NKERLSAFLCANADGTHKLSIIIGKSKLPK  
SVKEDTSTLPVIYKPSKDVWFTRELFSEWFFQNFVPE  
VRHFQLNVLRFHDEDVRA~~L~~LLDSCPAHPSESLS~~E~~D  
GRIKCMFFPHNTSTLIQPMNQGVILSCKRLYRWKQLE  
ESLVIFEESDDEQEKGDKGVSKIYNIKSAIFNWAKS  
WEEVKQITIANAWENLLYKKEPEYDFQGLEHGDYREIL  
EKCGELET~~K~~DDDRVWLN~~G~~DEEKGLLTKGGITKEV  
VQKGGAEKQTA~~E~~FKLSAVRESLDYLLDFV~~D~~ATPEFQ  
RFHFTLKEMQQEIVKQFQSKIHSRIGSFLKPRPHNIK  
DSFSGPSTSGSNH

(C-353aa encoded by BC069623)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.