FOR IN VITRO RESEARCH USE ONLY. NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ATRNL1 protein



Basic Information

Catalog Number:

Ag16270

Size:

50 µg Form:

Available lyophilized

Species: human

Expression Source:

e coli.-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

PYCKANCGSPDHGYCDLTGEKLCVCNDSWQGPDCS LNVPSTESYWILPNVKPFSPSVGRASHKAVLHGKFMW VIGGYTFNYSSFQMVLNYNLESSIWNVGTPSRGPLQR

(247-356 aa encoded by BC139916)

Selected Validation Data