

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# L-VOCC Fusion Protein



## Basic Information

### Catalog Number:

Ag16455

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Validated Application:

Blocking peptide

### Peptide Sequence:

EVKMNHDTEACSEPSLLSTEMLSYQDDENRQLTLPEE  
DKRDIRQSPKRGFLRSASLGRRASFLHLECLKRQKDRGG  
DISQKTVLPLHLVHHQALAVAGLSPLLQRSHSPASFPR  
PFATPPATPGSRGWPPQPVPTLRLEGVESSEKLNSSFP  
SIHCGSWAETTPGGGGSSAARRVRPVSLMVP SQAGA  
PGRQFHGSASSLVEAVLISEGLGQFAQDPKFIEVTTQ  
ELADACDMTIEEMESAADNLSGGAPQSPNGALLPFV  
NCRDAGQDRAGGEEDAGCVRARGRPSEELQDSRV  
YVSSL

(1835-2135 aa encoded by BC146846)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.