

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human HELT protein



## Basic Information

Catalog Number:

Ag16520

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

I L E M T V Q Y L R A L H S A D F P R G R E K E L L A E F A N Y F H Y G Y H  
E C M K N L V H Y L T T V E R M E T K D T K Y A R I L A F L Q S K A R L G A  
E P A F P P L G S L P E P D F S Y Q L H P A G P E F A G H S P G E A A V F  
P Q G S G A G P F W P P G A A R S P A L P Y L P S A P V P L A S P A Q  
Q H S P F L T P V Q G L D R H Y L N L I G H A H P N A L N L H T P Q H P P  
V L

(54-241 aa encoded by BC140795)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.