

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human IQCC protein



## Basic Information

**Catalog Number:**

Ag17030

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

CRDHSSWLQMKQNRKPSQEKRDTTRMENPEATDQ  
RLPHSQPQLQELQYHRSHLAMELLWLQQAINSRKEYL  
LLKQTLRSPEAGPIREEPRVFLEHGEQACERDQSQPS  
APLEDQSYDRRTTGELEQEDDSCHRVKSPHRSPGSLA  
TTQKNIA GAKCREPCYSKSGPPSSIPSNQALGDRLTK  
GPDDGRQTFGGTCLLQMKILEDQTPRGLKPRNHCPRK  
SRTQLSALYEDSNIKEMSPRKLHKPD CRTV RTQELG  
LSEDHIIWDGTLGGPEHSVLDLWRTKPPKGQAPTDRS  
SRDGTSNEPSHEGQKKQRTIPWRSKSPEILSSTKAGCT  
GEEQWRGRPWKTEPPG

(117-466 aa encoded by BC098317)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.