

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human VPS52 protein



Basic Information

Catalog Number:

Ag17611

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

EQRYPF~~EALFRSQHYALLDNSCREYLFICEFFV~~SGPA
AHDLFHAVMGR~~TL~~SMTLKHLDSYLADCYDAIAVFLCI
HIVLRF~~RNIAAKRDVPALDRYWEQVLALLWPRFELILE~~
MNVQSV~~RSTDPQRLGGLDTRPHYITRRYAEFSSALVSI~~
NQ~~TIPNERTMQLLGQLQVEVENFVLRVAAEFSSRKE~~
QLVFL~~INNYDMMGLVLMERAADDSKEVESFQQLLNA~~
RTQEF~~IEELLSPFGGLVAFVKEAEALIERGQAERLRGE~~
EARVT~~QLIRGFGSSWKSSVESLSQDVMRSFTNFRNGT~~
SIIQ~~GALTQLIQLYHRFHRVLSQPQLRALPARAELINIH~~
HLMVEL~~LKKHKPNF~~

(373-723 aa encoded by BC032108)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.