

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human SP1 protein



Basic Information

Catalog Number:

Ag17828

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6⁺His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSDQDHSMDEMTAVVKIEKGVGGNNGNGNGGGA
FSQARSSSTGSSSSSTGGGGQESQPSPLALLAATCSRI
ESPNENSNSQGPSQSGGTGELDLTATQLSQGANG
WQIISSSSGATPTSKEQSGSSTNGSNGSESSKNRTVS
GGQYVVAAPNLQNQVLTGLPGVMPNIQYQVIPQ
FQTVDGQQLQFAATGAQVQDGSQGIIPGANQ
QIITNRGSGGNIIAAMPNLLQQAAPLQGLANNVLSGQ
TQYVTNPVALNGNITLLPVNSVSAATLTPSSQAVTI
SSSGSQESGSQPVTSGMTISSASLVSSQASSSSFFTN
ANSYSTTTTNSMGMNFTTSGSSGNSQGGTPQRV
SGLQGS DALNIQNNQTS GGS LQAGQQKEGEQNNQ
TQQQIILIQPQLVQGGQALQALQAAPLSGQTFFTQ
AISQETLQNLQAVPNSGPIIIRTPTVGPNGQVSWQ
TLQLNLQVQNPQAQTITLAPMQGVSLGQTSSTNT
LTPIASAASIPAGTVTVNAAQLSSMPGLQTINLSALGT
SGIQVHPIQGLPLAIANAPGDHGAQLGLHGAGGDGI
HDDTAGGEEGENSPDAQPQAGRRTREACTCPYCK
DSEGRSGDPGKIKKHICHIQCGKIVYKTSHLRAHL
RWHTGERPFMCTWSYCGKRFTRDELQRHKRHTHTGE
KKFACPECPKRFMRSDHLSKHIKTHQNKKGPGVALS
VGTLPDLSGAGSESGTATPSALITTNMVAMEAICPE
GIARLANSGINVMQVADLQINISGNGF

(1-785aa encoded by BC062539)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.