

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PRDM14 protein



Basic Information

Catalog Number:

Ag18025

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6^{His}.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MALPRPSEAVPQDKVCYPPESSQNLAAYTFFPSYG
HYRNSLATVEEDFQPFQLEAAA SAAPAMPPFFRM
APLLSPGLGLQREPLYDLPWYSKLPWYPIPHVPREV
PPFLSSHEYAGASSED LGHQIIGGDNESG PCCGPD
LIPPPADASLLPEGLR TSQLLPSPSKQSE DGPKPSN
QEGKSPARFQFTEEDLHFVLYGVTPSLEHPASLHHAIS
GLLVPPDSSGSDSLPQTLDEDSLQ LPEGLCLMQTVFG
EVPHFGVFCSSFI AKGVRFGPFQ GKVVNASEVKTYG
DNSVMWEIFEDGHL SHFIDGKGGTGNWMSYVNCAR
FPKEQNLVAVQCQGHIFYESCKEIHQNOELLVWYGD
CYEKFLDIPVSLQVTEPGKQPSGPSEESAEGYRCERC
GKVFTYKYRDKHLKYTPCVDKGDRKFP CSLCKRSFEK
RDLRLRIHLHVHEKHRPHKSTCGKCFSSSLNKHMR
VHSGDRPYQCVYCTKRFTASSILRTHIRQHSGEKPFK
CKYCGKSFASHAAHDSHVRRSHKEDDGCSCSICGKIF
SDQETFYSHMKFHEDY

(1-571 aa encoded by BC052311)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.