

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human SORL1 protein



## Basic Information

Catalog Number:

Ag18277

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

PRNVVYGIFYATSFLDLYRNPISLTTSLHNKTVIVSKDE  
QYFLVLRVVVPYQGPSSDYVVVKMIPDSRLPPRHLHV  
VHTGKTSVVIKWESEPYDSDQDLYAIAVKDLIRKTD  
SYKVKSRNSTVEYTLNKLPGGKYHIIVQLGNMSKDSS  
IKITTVSLSAPDALIKIITENDHVLLFWKSLALKEKHF  
NESRGYEIHMFD SAMNITAYLGNTTDNFFKISNLKMG  
HNYTFTVQARCLFGNQICGEPAILLYDELGSGADASAT  
QARSTDVAAVVVPILFLLLSLGVGFALYTKHRRLLQ  
SSTAFANSHYSSRLGSAIFSSGDDLGEDEDEAPMITG  
FSDVPMVIA

(C-350aa encoded by BC137171)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.