

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# UPB1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag18996

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MAGAEWKSLEECLEKHLPLDLQEVKRVLYGKELRKLD  
LPREAFEAASREDFELQGYAFEAAEEQLRRPRIVHVGL  
VQNRIPANAPVAEQVSALHRRIKAIVEVAAMCGV  
NIICFQEAWTMPFAFCTREKLPWTEFAESAEDGPTTR  
FCQKLAKNHDMVVSPILERDSEHGDLWNTAVVIS  
NSGAVLGKTRKNHPRVGDFFNESTYYMEGNLGHVPV  
QTQFGRIAVNICYGRHPLNWLMSINGAEIIFNPSA  
TIGALSESLWPIEARNAAIANHCFTCAINRVGTEHFPN  
EFTSGDGKKAHQDFGYFYGSSVVAAPDSSRTPGLSRS  
RDGLLVAKLDLNLCCQQVNDVWNFKMTGRYEMYAREL  
AEAVKSNYSPTIVKE

(1-384aa encoded by BC131703)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.