

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human PAX2 protein



## Basic Information

**Catalog Number:**

Ag19008

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MDMHCKADPFAMHPGHGGVNQLGGVFVNGRPLP  
DVVRQRIVELAHQGVRPCDISRQLRVSHGCVSKILGR  
YYETGSIKPGVIGGSKPKVATPKVVDKIAEYKRQNPTM  
FAWEIRDLLAEGICDNDTVPSVSSINRIIRTKVQQPF  
HPTPDGAGTGVTAPGHTIVPSTASPPVSSASNDPVG  
SYSINGILGIPRSNGEKRRKRDVEVYTDPAHIRGGGGL  
HLVWTLRDVSEGSVPNGDSQSGVDSLKHLRADTFT  
QQQLEALDRVFERPSYDPVFQASEHIKSEQGNEYSLP  
ALTPGLDEVKSSLSASTNPELGSNVSGTQTPYVVVTGR  
DMASTTLPGYPPHV PPTGQGSYPTSTLAGMVPGSEFS  
GNPYSHQPQYTA YNEAWRF SNPALLMPPPGPLPLLP  
MTATSYRGDHIKIQADSFGLHIVPV  
(1-431aa encoded by BC148710)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.