FOR IN VITRO RESEARCH USE ONLY. NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human FAM20A protein



Basic Information

Catalog Number:

Ag19429

Size:

50 µg Form:

Available lyophilized

Species: human

Expression Source:

e coli.-derived, PGEX-4T, with N-terminal GST.

Biological Activity: Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

GFLIHLDNARGFGRHSHDEISILSPLSQCCMIKKKTLLH LQLLAQADYRLSDVMRESLLEDQLSPVLTEPHLLALDR RLQTILRTVEGCIVAHGQQSVIVDGPVEQSAPDSGQ

(424-541 aa encoded by BC136686)

Selected Validation Data