

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human PRDM4 protein



## Basic Information

**Catalog Number:**

Ag19675

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6<sup>°</sup>His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

EVAEWTDKAVNHIWIKIYHNGVLEFCIITDENEENW  
MMFVVRKARNREEQNLVAYPHDGKIFFCTSQDIPPENE  
LLFYYSRDYAAQIGVPEHPDVHLCNCGKECNSYTEFK  
AHLTSHIHNLPTQGHSGSHGSPHSKERKWKCSMCP  
QAFISPSKLVHFMGHMGMKPHKCDFCSKAFSDPSN  
LRTHLKIHTGQKNYRCTLCDKSFTQKAHLESHMVIHT  
GEKNLKC DYCDKLFMRQDLKQHVLIHTQERQIKCPK  
CDKFLRTNHLKHLNSHEGKR DYVCEKCTKAYLT KYH  
LTRHLKTCKGPTSSSSAPEEEEEEDDSEEDLADSVGTE  
DCRINSAVYSADELSA HK  
(451-801 aa encoded by BC035581)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.