

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human UHRF1BP1 protein



## Basic Information

Catalog Number:

Ag19683

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>+</sup>His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

LGRDRMTSTMHKMLSLPPAKEPMAKTDEGVAAPVSG  
GAARLRFFSMKRTVSQQSFDGVSLSGDPEDRISVDS  
DGSDSFVMLESESGPESVPPGSLSNVSDNAGVQGS  
PLVNNYGQGSPAANSSVSPGEDLIFHPVSVLVKVN  
EVSFGIEVRGDLTVALQAEELTLQLGTVGLWQFLH  
GQCPTGCFQESSTLKTGHIRPAVGLRFEVGPAAVH  
SPLASQNGFLHLLHGCDELLTSVLSGLPFLDEEIP  
VVVPMQIELNSSITLKDDIPPIYPTSPGPIPITLAMEHV  
VLKRSDDGVFHIGAAAQDKPSAEVLKSEKRQPPKEQ  
VFLVPTGEVFEQQT

(1053-1400aa encoded by BC137120)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.