

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# LGALS3BP Fusion Protein



## Basic Information

**Catalog Number:**

Ag19700

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6<sup>His</sup>.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

VNDGDMRLADGGATNQGRVEIFYRQWGTVCNDL  
WDLTDASVVCRALGFENATQALGRAAFQGSGPIML  
DEVQCTGTEASLADCKSLGWLKSNCRHERDAGVVCT  
NETRSTHTLDSRESEALGQIFDSQRGCDLSISVNVQ  
GEDALGFCGHTVILTANLEAQALWKEPGSNVTMSVD  
AECVPMVRDLLRYFYRRIDITLSSVKCFHKLASAYGA  
RQLQGYCASLFAILLPQDPSFQMPLDLYAYAVATGDA  
LLEKLCQLFLAWNFEALTQAEA WPSVPTDLLQLLPRS  
DLAVPSELALLKAVDTWSWGERASHEEVEGLVEKIRF  
PMMLPEELFELQFNLSLYWSHEALFQKKTLOALEFHT  
VPFQLLARYKGLNLTEDTYKPRIYTSPTWSAFVTDSSW  
SARKSQLVYQSRRGPLVKYSSDYFQAPSDYRYYPYQS  
FQTPQHPSFLFQDKRVSWSLVYLPTIQSCWNYGFSC  
SSDELVLGLTKSGGSDRTIAYENKALMLCEGLFVADV  
TDFEGWKAAPSAALDTNSSKSTSSFPCPAGHFNFRF  
VIRPFYLTNSSGVD

(19-585 aa encoded by BC002403)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.