

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human KAL1 protein



## Basic Information

**Catalog Number:**

Ag19718

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

FWSWMVSSKSLVPTKKIKRRKTTDGFQNSVILEKLQPD  
CDYVVELQAITYWGQTRLKSAKVS LHFTSTHATNNKE  
QLVKTRKGGIQTQLPFQRRRPTRPLEVGAFFYQDGQL  
QVKVYWKKTEDPTVNRVHVRWFPEACAHNRTTGSE  
ASSGMTHENYIILQDLSFSCYKVKTVQPIRPSKSHKAE  
AVFFTPPCSALKGKSHKPIGCLGEAGHVLSKVLAKPE  
NLSASFIVQDVNITGHFSWKMAKANLYQPMTGQVT  
WAEVTTESRQNSLPNSIISQSQILPSDHYVLTVPNLRP  
STLYRLEVQVLTGGEGPATIKTFRTPELPPSSAHRSHL  
KHRHPHHYKPSPERY

(331-680 aa encoded by BC137426)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.