

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human MAML2 protein



## Basic Information

**Catalog Number:**

Ag20163

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

VGNMQPTAQYSGGSSTISLNSNQALANPVSTHTILTP  
NSSLLSTSHGTRMPSTAVQNMGMVGNLPCNQPNT  
YSVTSGMNQLTQQRNPKQLLANQNNPMMRPPTLG  
PSNNNNVATFGAGSVGNSQQLRPNLTHSMASMPQ  
RTSNVMTSNTTAPNWASQEGTSKQQEALTSAGVRF  
PTGTPAA YTPNQLQAVGSQQF5QRAVAPPNQLT  
PAVQMRPMNQMSQTLNGQTMGPLRGLNLRPNQLST  
QILPNLNQSGTGLNQSRGTGINQPPSLTPSNFSPNQS  
SRAFQGTDHSSDLAFDFLSQQNDNMGPALNSDADFI  
DSLLKTEPGNDDWMKDINLDEILGNNS  
(805 - 1153 aa encoded by BC152449)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
[ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.