

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RECOMBINANT HUMAN AMH



## Basic Information

**Catalog Number:**  
Ag20166

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*E. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
MRDLPLTSLALVLSALGALLGTEALRAEEPAVGTSGLIFREDLDW  
PPGSQOEPLCLVALGGDSNGSSSPLRV/GALSAYEQAF LGAVQ  
RARWGPRDLATFGV CNTGDRQAALPSLRRLGAWLRDPGGQR  
LVLHLEEVTWEPTPSLRFQEPFGGAGPPELALLVLYPGGPE  
VTVTRAGLPGAQSLCPSRDTRYLVAVDRPAGAWRGSGLALTL  
QPRGEDSR LSTARLQALLFGDDHRCFTRMTPALLLPRSEPAFL  
PAHQQLDVTVPFPPRPSAELEESPPSADPFLETLTRLVRLRVPP  
ARASAPRLALDPDALAGFPQGLVNLSDPAALERLLDGEEPLLLLL  
RPTAATTGDPAPLHDPTSAPWATALARRVAAELQAAAAELRSLP  
GLPPATAPLLARLLALCPGGPGGLGDLRALLL KALQGLRVEW  
RCRDRPGPGRAGRSAGATAADGPCALRELSVDLRAERSVLIFE  
TYQANNCCQVCGWPGSDRNPRYGNHVLLLLKMQARGAALAR  
PPCCVPTAYAGKLLLSLSEERISAHMVPNMWATECGCR  
(1-560 aa encoded by BC049194)

## Selected Validation Data

For technical support and original validation data for this product, please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),  
or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.