

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Ubiquilin 2 Fusion Protein



## Basic Information

**Catalog Number:**

Ag20251

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6<sup>His</sup>.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MAENGESSGPPRPSRGPAAQGSAAAPAEPKIIKVTV  
KTPKEKEEFAVPENSSVQQFKEAISKRFKSQTDQLVLI  
FAGKILKDQDTLIQHGIHDGLTVHLVIKSQNRPGQS  
TQPSNAAAGNTTTSASTPRSNSTPISTNSNPFGLGSLG  
GLAGLSSLGLSSTNFSELQSQQMQLMASPEMMIQI  
MENPFVQSMLSNPDLMRQLIMANPQMQQLIQRNPEI  
SHLLNPDIMRQTLIARNPAMMQEMMRNQDLALSN  
LESIPGGYNALRRMYTDIQEPLNAAQEQFGGNPFAS  
VGSSSSSGEGTQPSRTENRDPLPNPWAPPATQSSA  
TTSTTTSTGSGSGNSSNATGNTVAAAANYVASIFSTP  
GMQSLQQTENPQLIQNMLSAPYMRSMQSLSQNP  
DLAAQMMLNSPLFTANPQLQEQMRPQLPAFLQMQ  
NPDTLSAMSNPRAMQALMQIQQGLQTLATEAPGLIP  
SFTPGVGVVLGTAIGPVGPVTPIGPIVFPFTPIGPI  
GPIGPTGPAAPPSTGSGGPTGPTVSSAAPSETTSPT  
SESGPNQFIIQQMVQALAGANAPQLPNPEVRFQQ  
QLEQLNAMGFLNREANLQALIATGGDINAAIERLLGS  
QPS

(1-624 aa encoded by BC069237)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.