

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human JNK protein



## Basic Information

Catalog Number:

Ag21426

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET30a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSRSKRDNNFYSVEIGDSTFTVLKRYQNLKPIGSGAQ  
GIVCAAYDAILERNVAIKKLSRPFQNTAKRAYRELV  
LMKCVNHKNII GLLNVFTPQKSLEEFQDVYIVMELMD  
ANLCQVIQMELDHERMSYLLYQMLCGIKHLHSAGIIH  
RDLKPSNIVVKS DCTLKLDFGLARTAGTSFMMTPYVV  
TRYRAPEVILGMGYKENVDIWSVGCIMGEMIKGGVL  
FPGTDHIDQWNKVIEQLGTPCPEFMKKLQPTVRTYVE  
NRPKYAGYSFEKLPDVLFPADSEHNKLKASQARDLLS  
KMLVIDASKRISVDEALQHPYINVWYDPSEAEAPPPKI  
PDKQLDEREHTIEEWKELIYKEVMDLEERTKNGVIRGQ  
PSPLAQVQQ

(1-384 aa encoded by BC130572)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.