

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# PTPN11 Fusion Protein



## Basic Information

**Catalog Number:**

Ag21472

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MTSRRWFHPNITGVEAENLLTRGVDGSFLARPSKSN  
PGDFTLSVRRNGAVTHIKIQNTGDYDLYGGEKFATL  
AELVQYYMEHHGQLKEKNGDVIELKYPLNCADPTSER  
WFHGLSGKEAEKLLTEKKGHSFLVRESQSHPGDFV  
LSVRTGDDKGESNDGKSKVTHVMIRCQELKYDVGGG  
ERFDSLTDLVEHYKKNPMVETLGTVLQKQPLNTRIN  
AAEIESRVRELSKLAETTDKVKQGFWEFEFETLQQQEC  
KLLYSRIKEGQRQENKNKNRYKNI LPFDHTRVVLHDGD  
PNEPVS DYINANIIMPEFETKCNNSKPKSYIATQGCL  
QNTVNDVWRMVFQENS RVIVMTTKEVERGKSKCVK  
YWPDEYALKEYGVMRVRNVKESAAHDYTLRELKLSKV  
GQGN TERTVWQYHFRTWPDHGVPSDPGGVLDLFLEE  
VHHKQESIMDAGPVV VHC R  
(1-460 aa encoded by BC008692)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.