

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

LCK Fusion Protein



Basic Information

Catalog Number:

Ag21475

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MGCGCSSHPEDDWMENIDVCENCHYPIVPLDGKGTLLIRNGSEVRDPLVYEGSNPPASPLQDNLVIALHSYEP SHDGD LGFEKGEQLRILEQSGEWWKAQSLTTGQEGF IPFNFAKANSLEPEPWFFKNLSRKDAERQLLAPGNT HGSFLIRESESTAGSFSLSVRDFDQNGGEVVKHYKIR NLDNNGGFYISPRITFPGLHELVRHYTNASDGLCTRLSRP CQTQKPQKPWWEDEWEVPRET LKLVRLGAGQFGE VWMGYNGHTKVAVKSLKQGSMSPD AFLAEANLMK QLQHQRRLRYAVVTQEPIYIITEYMENDTLLDSQLEE KGLGASPWNLGQQLLLPTGSLVDFLKTPSGIKLTINK LLDMAAQIAEGMAFIEERNYIHRDLRAANILVSDTLSC KIADFGLARLIEDNEYTAREGAKFKIKWTAPEAINVGT FTIKSDVWSFGILLTEIVTHGRIPYPGMTNPEVIQNLER GYRMVRPDNCPPELYQLMRLCWKERPEDRPTFDYLRSLVLEDDFTATEGQYQPQP

(1-539 aa encoded by BC013200)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.