

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# FAM105A Fusion Protein



## Basic Information

**Catalog Number:**

Ag22243

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PET30a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MAATRSPTRARERERSGAPAGSDQVHSWMLATSQ  
ALDTVWRMAKGFVMLAVSFLVAAICYFRRLHLYSGHK  
LKWWIGYLQRKFKRNLSVEAEVDLLSYCAREWKGETP  
RNKLMRKA YEELFWRHHIKCVRQVRRDNYDALRSVLF  
QIFSQGISFSPSWMKEKDIVKLPEKLLFSQGCNWIQQY  
SFGPEKYTGSNVFGKLRKYVELLKTQWTEFNIGIRDYHK  
RGS MCNTLFSDAILEYKLYEALKFIMLYQVTEVYEQMK  
TKKVI PSLFRLFSRETSSDPLSFMMNHLNSVGDTCGL  
EQIDMFILGYSLEVKIKVFRLFKFNRSRDFEVCYPEEPLR  
DWPEISLLTENDRHYHIPVF

(1-356 aa encoded by BC011524)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.