

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human BRUNOL4 protein



## Basic Information

Catalog Number:

Ag23043

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

LKPLFEFGIKIYELTVLKDRFTGMHKGCAFLTYCERESA  
LKAQSALHEQKTLPGMNRPIQVKPADSESRRGGSSCLR  
QPPSHRKLFGMLNKQSQSEDDVRRLEAFGNIECTI  
LRGPDGNSKGCFAVKYSSHAEQAALHGSQTMP  
GASSSLVVKFADTKERTMRRMQQMAGQMGMFNP  
MAIPFGAYGAYQALMQQAALMASVAQGGYLN  
MAFAAAQMQMAALNMNGLAAAPMTPTS  
(70-315 aa encoded by BC004167)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.