

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human P85a; PIK3R1 protein



Basic Information

Catalog Number:

Ag2344

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

RDTADGTFLLVRDASTKMHGDTLTLRKGGNNKLIKIFH
RDGKYGFSDDLTFSSVVELINHYRNESLAQYNPKLDVK
LLYPVSKYQQDQVVKEDNIEAVGKKLHKYNTQFQEK
SREYDRLYEYTRTSQEIQMKRTAIEAFNETIKIFEEQC
QTQERYSKYIEKFKREGNEKEIQRMHNYDKLSRISE
IIDSRRRLEEDLKKQAAEYREIDKRMNSIKPDLIQLRKTR
DQYLMWLTQKGVQRQKKLNEWLGNTEDQYSLVED
DEDLPHHDEKTWNVGSSNRNKAENLLRGKRDGTFLLVR
ESSKQGCYACSVVVDGEVKHCVINKTATGYGFAEPY
NLYSSLKELVLHYQHTSLVQHNDLSLNTLAYPVYAQQ
RR

(78-454 aa encoded by BC030815)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.